

RESEARCH ARTICLE OPEN ACCESS

CRISPR-Enhanced RAA-SHERLOCK Assay for Point-of-Care Detection of Cyprinid Herpesvirus-3: Development, Validation and Clinical Application

Hatem Soliman^{1,2,3} | Naveed Akram¹ | Mona Saleh¹¹Clinical Division of Fish Health, University of Veterinary Medicine, Vienna, Austria | ²Academy of Scientific Research & Technology, Cairo, Egypt | ³Department of Aquatic Animal Medicine and Management, Faculty of Veterinary Medicine, Assiut University, Assiut, Egypt**Correspondence:** Mona Saleh (mona.saleh@vetmeduni.ac.at)**Received:** 27 August 2025 | **Revised:** 15 October 2025 | **Accepted:** 18 October 2025**Funding:** Hatem Soliman gratefully acknowledges the financial support provided by the Academy of Scientific Research and Technology (ASRT), Egypt, through the ASRT-STAR fellowship program, Cycle 1. This work was supported in part by the Veterinary Medicine University, Vienna.**Keywords:** CRISPR-Cas13a | CyHV-3 | lateral flow detection | point-of-care diagnostics | recombinase-aided amplification | SHERLOCK

ABSTRACT

Koi herpesvirus disease (KHVD), caused by Cyprinid herpesvirus-3 (CyHV-3), poses a significant threat to global aquaculture due to its high mortality rates and economic impact. Current diagnostic methods, such as PCR, are limited by equipment dependency and procedural complexity, hindering point-of-care (POC) applications. To address this, we developed an integrated assay combining recombinase-aided amplification (RAA) with CRISPR-Cas13a-mediated SHERLOCK technology and lateral flow detection (LFD) for rapid and visual detection of CyHV-3 in clinical samples. The KHV-SHERLOCK assay targets a conserved region of the CyHV-3 thymidine kinase (TK) gene, demonstrating exceptional specificity with no cross-reactivity to related pathogens or host DNA. Sensitivity evaluations revealed a detection limit of 100 ag/μL for CyHV-3 plasmid DNA, tenfold more sensitive than the conventional PCR (1 fg/μL) assay, even in the presence of 100 ng of carp genomic DNA as background interference. Clinical validation using 50 archived samples showed 100% concordance with reference PCR results, confirming diagnostic reliability. The assay's isothermal RAA step (37°C, 40 min) and CRISPR-Cas13a detection (37°C, 1 h) enable equipment-free operation, while LFD provides unambiguous visual results within minutes. This platform merges high sensitivity with POC practicality, offering a transformative tool for field-based KHVD surveillance.

1 | Introduction

Koi herpesvirus disease (KHVD), caused by Cyprinid herpesvirus 3 (CyHV-3; initially classified by Waltzek et al. (2005) and reclassified as Cyvirus cyprinidallo 3 in the ICTV 2022 taxonomy release), is a highly contagious viral infection affecting both common carp (*Cyprinus carpio carpio*) and ornamental koi (*Cyprinus carpio koi*) (Hedrick et al. 2000; Aoki et al. 2007). The disease has caused significant economic losses in global aquaculture due to its high virulence and capacity to induce latent infections (Ronen et al. 2003; Haenen et al. 2004;

Gotesman et al. 2013). Clinical manifestations of KHVD include severe gill necrosis and systemic pathology, with mortality rates reaching 70%–100% within the permissive temperature range of 18°C–28°C (Hedrick et al. 2000; Gilad et al. 2003; Perelberg et al. 2003).

Although virus isolation is considered the gold standard for viral infections diagnostics, its application for KHVD is limited by the restricted availability of permissive cell lines and lower sensitivity compared to molecular methods (Haenen et al. 2004; World Organization for Animal

This is an open access article under the terms of the [Creative Commons Attribution](https://creativecommons.org/licenses/by/4.0/) License, which permits use, distribution and reproduction in any medium, provided the original work is properly cited.

© 2025 The Author(s). *Journal of Fish Diseases* published by John Wiley & Sons Ltd.

Health [WOAH] (2022)). Consequently, PCR-based techniques, including real-time PCR (Gilad et al. 2004) and conventional PCR (Bercovier et al. 2005), are the recommended approaches for presumptive diagnosis in clinically affected fish. Additional PCR assays (Yuasa et al. 2005; El-Matbouli et al. 2007; Bergmann et al. 2010; Engelsma et al. 2013), isothermal nucleic acid amplification methods (Gunimaladevi et al. 2004; Yoshino et al. 2009; Soliman and El-Matbouli 2005, 2010, 2018; Prescott et al. 2016; Li, Zhang, et al. 2024; Li, Li, et al. 2024) and gold nanoprobe-based assays (Saleh and El-Matbouli 2015; Bavarsad et al. 2025), offer comparable or superior sensitivity for CyHV-3 detection.

Recombinase aided amplification (RAA), an advanced derivative of recombinase polymerase amplification (RPA), has been widely adopted for sensitive detection of bacterial, viral and parasitic pathogens (Wang et al. 2022; Zhao et al. 2023; Yan et al. 2023). Compared to conventional PCR, RAA demonstrates significant practical advantages, such as reduced reliance on expensive thermal cycling equipment and specialised consumables (Li et al. 2023). This isothermal amplification method offers broader application potential while maintaining high detection accuracy (Bei et al. 2010). RAA systems exhibit superior energy efficiency and operational simplicity, making them particularly suitable for field applications (Wang et al. 2022). Recent advances have established RAA as an emerging nucleic acid amplification platform for various aquatic pathogen detections with validated applications (Gui et al. 2022).

Recently, CRISPR-Cas systems have emerged as a transformative tool for nucleic acid detection, offering rapid, highly specific and efficient diagnostic capabilities, thereby revolutionising the field of molecular diagnostics (Hatoum-Aslan 2018; Yin et al. 2021). Several CRISPR-Cas nucleases, such as Cas12a, Cas12b and Cas13a, possess dual functionality: they cleave target-specific nucleic acids and exhibit collateral cleavage activity, resulting in nonspecific degradation of nearby single-stranded DNA or RNA upon target recognition (Abudayyeh et al. 2016; Chen et al. 2018; Li et al. 2018, 2019). This property has been exploited to develop highly sensitive and specific diagnostic platforms for viral detection (Gootenberg et al. 2017; Chen et al. 2018). A prominent example is SHERLOCK (Specific High-Sensitivity Enzymatic Reporter UnLOCKing), which combines Cas13a's collateral cleavage activity with isothermal amplification for robust nucleic acid detection (Gootenberg et al. 2017; Kellner et al. 2019).

The SHERLOCK platform serves primarily as a qualitative diagnostic system, providing a binary readout of target nucleic acid presence or absence (Gootenberg et al. 2017). While demonstrating exceptional sensitivity with detection limits reaching attomolar concentrations (Chen et al. 2018), the platform lacks the quantitative capacity of real-time PCR for precise viral load determination (Kellner et al. 2019). Recent methodological advancements incorporating fluorescence or lateral flow detection (LFD) have enabled semi-quantitative analysis in optimised systems (Myhrvold et al. 2018). The platform has demonstrated broad diagnostic utility, successfully detecting diverse pathogens including White Spot Syndrome Virus, Dengue virus, Zika virus, Infectious

Bursal Disease Virus, *Yersinia ruckeri*, *Neisseria gonorrhoeae* and *Plasmodium* species (Myhrvold et al. 2018; Sullivan et al. 2019; Lee et al. 2020; Allan-Blitz et al. 2023; Calderón et al. 2024; Ramadan et al. 2025). SHERLOCK's modular crRNA design enables rapid system reconfiguration, offering particular advantages for emerging pathogen surveillance and strain-specific differentiation in outbreak settings (Ramadan et al. 2025).

In this study, we developed a field-deployable diagnostic assay by integrating RAA with the SHERLOCK platform and LFD for rapid, sensitive and specific visual detection of CyHV-3 DNA in fish tissues. We systematically evaluated the assay's specificity, sensitivity and clinical applicability, highlighting its potential for point-of-care (POC) diagnosis.

2 | Materials and Methods

2.1 | Nucleic Acid Templates

Genomic DNA was extracted from purified CyHV-3 using the QIAamp DNA Mini Kit (QIAGEN GmbH, Germany), according to the manufacturer's instructions. To evaluate the specificity of the assay, genomic DNA from specific pathogen-free (SPF) common carp (ZODIAC, Wageningen, the Netherlands) and non-infected common carp brain (CCB) cells was included as negative controls. Additionally, genomic DNA from Cyprinid herpesvirus-1 (CyHV-1) and from Cyprinid herpesvirus-2 (CyHV-2), generously provided by Dr. Bergmann (Friedrich-Loeffler-Institute, Greifswald-Insel Riems, Germany), was used to assess cross-reactivity. Further specificity testing was conducted using cDNA from spring viraemia of carp virus (SVCV) and DNA from a confirmed carp edema virus (CEV)-positive sample to validate the discriminatory capacity of the newly developed KHV-SHERLOCK assay.

2.2 | Cloning and Sequencing of a CyHV-3 Thymidine Kinase (TK) Gene Fragment

A 409-bp fragment of the CyHV-3 thymidine kinase (TK) gene was amplified by PCR using 250 ng of genomic DNA and the primer pair KHV-TKI/KHV-TKr, as previously described by Bercovier et al. (2005). The resulting PCR products were ligated into the pCR4-TOPO vector (Thermo Fisher Scientific) and transformed into One Shot TOP10 chemically competent *Escherichia coli* using the TOPO TA Cloning Kit for sequencing (Thermo Fisher Scientific), in accordance with the manufacturer's protocol.

Recombinant plasmids were isolated using the QIAprep Spin Miniprep Kit (QIAGEN GmbH) following the manufacturer's instructions. The inserted TK gene fragments were subsequently sequenced (Eurofins Genomics, Germany), and the obtained sequences were aligned and compared with the reference sequence from Bercovier et al. (2005) to confirm accuracy. For downstream applications, the plasmids were linearised with the *NotI* restriction enzyme (Thermo Fisher Scientific, Germany) and purified using the QIAquick PCR Purification Kit (QIAGEN GmbH).

2.3 | Design of KHV-RAA Primers, -crRNA Sequences and Reporter Probe

The nucleotide sequences of the CyHV-3-TK gene available in the GenBank were retrieved and aligned using ClustalW (Thompson et al. 1994) to identify conserved regions. KHV-RAA primers (Table 1) were designed to target a conserved 139-nucleotide region of the CyHV-3-TK gene (GenBank accession: [AJ535112](#)), following the design criteria specified in the TwistAmp Assay Design Manual. Primer specificity was verified through in silico analysis using BLAST against the GenBank database. To facilitate downstream SHERLOCK detection, a T7 promoter sequence was incorporated at the 5' terminus of the forward KHV-RAA primer, enabling in vitro transcription by T7 RNA polymerase.

The KHV-specific crRNA was designed to comprise two functional elements: a 23-nucleotide spacer sequence targeting the amplified RAA product region between the primer binding sites and the *Leptotrichia wadei* Cas13a (LwaCas13a) direct repeat sequence. The spacer sequence was carefully selected to avoid complementarity with the KHV-RAA primer sequences, thereby minimising potential false-positive signals. Specificity of the spacer sequence for CyHV-3 was confirmed through comprehensive BLAST analysis. The final crRNA construct was generated by appending the LwaCas13a direct repeat sequence to the 5' end of the reverse complement of the selected spacer sequence (Table 1).

For visual detection of the SCHERLOCK assay results, a 14-nucleotide polyuridine (poly-U) single-stranded RNA LF-reporter probe was employed. The probe was dually labelled with 6-carboxyfluorescein (FAM) at the 5' end and biotin at the 3' end (Table 1), enabling LFD readout of the KHV-SHERLOCK assay results.

2.4 | Establishment of KHV-RAA Assay

The KHV-RAA assay was conducted using the RAA Nucleic Acid Amplification Kit (Xpedite Diagnostics, Hallbergmoos,

Germany) with minor modifications to the manufacturer's protocol. The reaction mixture was prepared on ice with the following composition: 25 µL of rehydration buffer, 2.5 µL of each forward and reverse KHV-RAA primer (10 µM each), and 13 µL of nuclease-free water. This mixture was then used to rehydrate the lyophilised RAA pellet. Subsequently, 2 µL of DNA template and 5 µL of magnesium acetate were added to the reaction mixture, followed by gentle pipetting to ensure homogenisation. The reaction was initiated by incubation at 37°C for 4 min. After brief manual mixing by inversion and centrifugation, the tubes were returned to the thermal block for a 40-min amplification period. A no-template control (NTC) was included in each experimental run to monitor potential contamination. The resulting KHV-RAA amplicons were either immediately processed for LwaCas13a-based detection or stored at -20°C for subsequent analysis.

2.5 | Establishment of the KHV-SCHERLOCK Assay

The 139bp RAA-amplified fragment of the CyHV-3-TK gene was transcribed into RNA using T7 RNA polymerase, and the resulting target RNA was subsequently detected via LwaCas13a-mediated collateral RNA cleavage, which induces the release of a fluorescent reporter signal. Detection was performed using an endpoint assay in which lateral flow strips were exposed to the reaction mixture following incubation. The Cas13a detection reaction was carried out on ice in a total volume of 20 µL, consisting of 11.12 µL of nuclease-free water (Thermo Fisher Scientific, Austria), 0.4 µL of HEPES (1M, pH 6.8), 0.18 µL of MgCl₂ (1M; Thermo Fisher Scientific, Austria), 2 µL of LwaCas13a nuclease (63.3 µg/mL; GenScript Biotech, Rijswijk, the Netherlands), 0.8 µL of rNTP mixture (25 mM; New England BioLabs, Germany), 1 µL of murine RNase inhibitor (40 U/µL; New England BioLabs, Germany), 0.5 µL of T7 RNA polymerase (20 U/µL; Promega Biotech, Denmark), 1 µL of KHV-crRNA (10 ng/µL), 2 µL of LF-RNA reporter probe (10 µM) and 1 µL of KHV-RAA amplicon. After thorough mixing, the reaction mixture was incubated at 37°C for 1 h and then immediately placed on ice to terminate the reaction prior to signal readout.

TABLE 1 | Oligonucleotide sequences used in this study, designed based on the CyHV-3 thymidine kinase (TK) gene sequence (GenBank accession no. [AJ535112](#)).

Sequence name	Oligonucleotides sequences (5'→3')	References
KHV-RAA-for ^a	5'- <i>GAAATTAATACGACTCACTATAGGG</i> GGACTTTATGCAGCAGCCCTTCAAGCAGGTGACG-3'	This study
KHV-RAA-Rev	5'- GTCCGTGCCCTGAGAGATTCTGACGGTGAAG-3'	This study
KHV-crRNA ^b	5'- <i>GAUUUAGACUACCCCAAAAACGAAGGGGACU</i> <i>AAAAC</i> AGCUUGUCCAGCUUGUCCGCCAU-3'	This study
LF-RNA-reporter ^c	5'- (6-FAM) – rUrUrUrUrUrUrUrUrUrUrUrU–BioTEG - 3'	Kellner et al. (2019)
KHV-TK-For	5'- GGGTTACCTGTACGAG -3'	Bercovier et al. (2005)
KHV-TK-Rev	5'- CACCCAGTAGAT TATGC -3'	

^aThe sequence in bold italics is the T7 RNA polymerase promoter sequence.

^bThe sequence in bold italics is the LwaCas13a direct repeat sequence.

^c6-FAM, 6-carboxyfluorescein; Bio TEG, biotin triethylene glycol spacer.

2.6 | Interpretation of KHV-SHERLOCK Assay Results

The KHV-SHERLOCK assay results were analysed using the HybriDetect 1–Universal Lateral Flow Assay Kit (Milenia Biotec GmbH, Germany). For detection, 20 µL of the KHV-SHERLOCK reaction products were mixed with 80 µL of HybriDetect assay buffer and applied to lateral flow strips. Results were visually assessed within 5 min of application.

In negative samples lacking the target gene, trans-cleavage activity does not occur, leaving the LF-reporter intact. In this case, the biotinylated end of the reporter is captured by streptavidin at the control line, while anti-FAM antibody-conjugated gold nanoparticles bind to the intact FAM moiety, producing a distinct purple control band. The absence of a test line under these conditions confirms negative results, indicating unsuccessful target gene amplification.

Conversely, in positive samples containing the target CyHV-3-TK gene, Cas13a-mediated trans-cleavage of the LF-reporter occurs. This cleavage results in reporter fragments where biotin, rather than FAM, becomes the predominant terminal moiety. Consequently, fewer intact FAM-labelled reporters reach the control line, while cleaved biotinylated fragments are captured at the test line by streptavidin. This produces two possible positive outcomes: a single, strong test line (indicating high target concentration) or both test and control lines (indicating moderate target concentration). The appearance of either pattern confirms the presence of the CyHV-3-TK target gene in the tested sample.

2.7 | Evaluation of the KHV-SCHERLOCK Assay Specificity

To evaluate the analytical specificity of the KHV-SCHERLOCK assay developed in this study, nucleic acids from CyHV-1, CyHV-2, SVCV and CEV were tested as templates under optimised assay conditions. To control for non-specific amplification, genomic DNA from SPF common carp and non-infected CCB cells was included, along with a NTC. All experiments were conducted in triplicate to ensure reproducibility and reliability of the results.

2.8 | Evaluation of the KHV-SCHERLOCK Assay Sensitivity

To evaluate the sensitivity of the KHV-SCHERLOCK assay, 10-fold serial dilutions of a KHV-TK gene plasmid (1 ng/µL; 2.1×10^8 copies/µL) were prepared and tested both in isolation and after pre-mixing with 100 ng of SPF common carp genomic DNA. Each reaction included a NTC to monitor potential contamination or non-specific amplification. The resulting KHV-SCHERLOCK amplicons were immediately assessed using lateral flow dipstick (LFD) analysis. The lower detection limit of the KHV-SCHERLOCK assay was compared against the conventional diagnostic PCR assay (Bercovier et al. 2005). All experimental procedures were conducted under previously established optimal assay conditions. To confirm reproducibility, the entire assay was independently replicated three times.

2.9 | Application of the KHV-SCHERLOCK Assay in Clinical Samples

To validate the diagnostic utility of the KHV-SCHERLOCK assay, the previously optimised protocol was applied to detect CyHV-3 in clinical specimens. A total of 50 archived fish DNA samples, previously characterised using conventional diagnostic methods, were analysed using the KHV-SCHERLOCK assay. The resulting detection outcomes were systematically compared with those obtained from standard reference assays to assess concordance and diagnostic performance.

3 | Results

3.1 | Optimisation and Performance of the KHV-SHERLOCK Diagnostic Assay

The KHV-SHERLOCK assay was developed as a two-step diagnostic procedure, initiated by the RAA assay followed by CRISPR-Cas13a-mediated detection. The KHV-RAA step was optimised using primers specifically targeting the CyHV-3-TK gene to ensure high specificity. Standard KHV-RAA reactions were performed at 37°C for 40 min in a heating block, with successful amplification of the target CyHV-3-DNA confirmed by agarose gel electrophoresis of purified amplicons. No amplification was observed in the NTC, confirming the absence of non-specific amplification (data not shown).

Following RAA, the SHERLOCK reaction mixture was applied to the amplification products, with optimal detection conditions determined to be 37°C for 1 h. The assay results were visualised using a LFD, where the positive control exhibited a clear test line band, confirming target detection, while the NTC displayed only a control line band, demonstrating assay specificity (Figure 1).

3.2 | Specificity of the Established KHV-SHERLOCK Assay

The analytical specificity of the KHV-SHERLOCK assay was rigorously assessed. The assay exhibited high specificity, as evidenced by a distinct single band on the LFD test line exclusively in the CyHV-3-positive control sample (Figure 2), confirming a true-positive detection. In contrast, no cross-reactivity was observed with genomic DNA extracted from SPF common carp, non-infected CCB cells or NTCs. Furthermore, the assay demonstrated no amplification or detection of closely related viral pathogens, including CyHV-1, CyHV-2, SVCV or CEV, as indicated by the absence of test line bands and the presence of only the control line band in these samples. These findings confirm the specificity of the KHV-SHERLOCK assay, ensuring accurate and exclusive identification of CyHV-3 without cross-detection of non-targeted pathogens.

3.3 | Sensitivity of the KHV-SHERLOCK Assay

The analytical sensitivity of the KHV-SHERLOCK assay was systematically compared to conventional diagnostic PCR using




FIGURE 1 | Visual detection of the KHV-SHERLOCK assay using lateral flow dipsticks. The lateral flow dipstick (LFD) displays distinct colorimetric signals corresponding to the assay results. C-line=The control line—the appearance of a dark purple band indicates a negative result for the target sequence. T-line=The test line—exhibits a dark purple signal only in the presence of the target CyHV-3- DNA, denoting a positive result. NTC=no-template control. +VECO=positive control.

10-fold serial dilutions of a KHV-TK gene plasmid (initial concentration: 1 ng/ μ L; 2.1×10^8 copies). To evaluate clinical relevance, all assays were performed both with and without 100 ng of SPF-common carp genomic DNA as a background matrix. Under optimised reaction conditions, the KHV-SHERLOCK assay demonstrated consistently detection down to 21 plasmid copies (100 ag/ μ L), regardless of background DNA presence. This represents a 10-fold improvement over conventional PCR, which exhibited a detection limit of 210 copies (1 fg/ μ L) under identical conditions (Figure 3).

3.4 | Evaluation of the KHV-SHERLOCK Assay Using Clinical Specimens

To validate the diagnostic utility of the KHV-SHERLOCK assay for CyHV-3 detection in clinical settings, a panel of 50 archived fish DNA samples, previously characterised using standard reference PCR methods, was analysed. The assay demonstrated complete concordance with the established diagnostic results, correctly identifying all CyHV-3-positive samples without any false-positive or false-negative detections. These findings confirm that the KHV-SHERLOCK assay exhibits high diagnostic accuracy and reliability when applied to clinical specimens.




FIGURE 2 | Specificity assessment of the KHV-SHERLOCK assay. The KHV-SHERLOCK assay specifically detects KHV-DNA, as indicated by a dark purple test line on lateral flow dipsticks. A control line signal indicates a negative result (no CyHV-3-DNA). No cross-reactivity was observed with Cyprinid Herpesvirus-1 (CyHV-1), Cyprinid herpesvirus-2 (CyHV-2), Cyprinid herpesvirus-3 (CyHV-3), Carp edema virus (CEV) or Spring viraemia of carp virus (SVCV). Genomic DNA from specific pathogen-free common carp (Carp DNA) and non-infected common carp brain cells (CCB DNA) yielded negative results. No-template controls (NTC) confirmed absence of contamination.

4 | Discussion

So far, no effective therapeutic agents have been developed to treat CyHV-3 infections in fish; accuracy and rapid detection of CyHV-3-DNA in host tissues are essential for disease control and prevention. The ability to detect CyHV-3 with high sensitivity and specificity is particularly critical in the carp and koi




FIGURE 3 | Sensitivity of the KHV-SHERLOCK assay. The assay demonstrates superior detection sensitivity compared to conventional diagnostic PCR. A 10-fold serial dilution series (lanes 1–9) of CyHV-3-plasmid DNA (initial concentration: 1 ng; 2.1×10^8 copies) was analysed. The KHV-SHERLOCK assay achieved a detection limit of 21 CyHV-3-plasmid copies (100 ag/ μ L), visualised by distinct purple bands at both the test (T-line) and control (C-line) positions on lateral flow dipsticks. In contrast, the diagnostic PCR assay detected only 210 CyHV-3-plasmid copies (1 fg/ μ L). Negative results were confirmed by a solitary control line signal. No-template controls (NTC) verified assay specificity, and a 100-bp DNA ladder (M) was included for reference.

aquaculture industries, where early diagnosis can mitigate large-scale outbreaks and minimise economic losses.

Currently, several nucleic acid-based detection techniques have been established for CyHV-3 identification, including conventional PCR and real-time PCR (Gilad et al. 2004; Bercovier et al. 2005). Although these techniques exhibit high analytical sensitivity, these methods are time-consuming, require expensive equipment, and are therefore unsuitable for small-scale clinical settings or narrow-range disease surveillance programmes. Their procedural complexity limits their practicality for POC applications, where rapid, user-friendly and cost-effective diagnostics are essential.

To overcome these limitations, we employed SHERLOCK, a next-generation nucleic acid detection technology, to develop an ultra-sensitive, field-deployable assay for CyHV-3 DNA detection in fish tissue. This approach combines high diagnostic sensitivity with the simplicity required for POC testing (Kellner et al. 2019), offering a promising solution for rapid and accessible KHV surveillance.

The KHV-SHERLOCK assay employs an isothermal preamplification step for CyHV-3-DNA using the RAA assay, followed by transcription of the amplified products via T7 RNA polymerase and subsequent detection using LwaCas13a with a lateral flow colorimetric readout. RAA was selected for preamplification due to its high sensitivity, specificity, ability to function at a constant low temperature, and rapid detection of low nucleic acid concentrations (Wang et al. 2019). Unlike conventional amplification methods, the RAA assay operates isothermally, completing amplification within 10–40 min without requiring stringent laboratory conditions, specialised expertise or costly instrumentation (Fan et al. 2019). Furthermore, its compatibility with portable devices, owing to the absence of thermal cycling, enhances its applicability in field settings (Yan et al. 2023).

In this study, the CyHV-3-TK gene was selected as the target for designing KHV-RAA primers and crRNA due to its critical role in CyHV-3 virulence, a feature conserved among herpesviruses (Boivin et al. 2002; Bercovier et al. 2005). Additionally, this gene is routinely targeted in established KHV-diagnostic PCR assays, allowing for direct sensitivity comparisons between conventional PCR and the SHERLOCK assay.

The RAA primers were designed to target a conserved region of the CyHV-3-TK gene, ensuring broad detection capability across all known CyHV-3 isolates. To maximise amplification efficiency, the primers were constructed with the minimal length required for recombinase-mediated integration into duplex DNA, facilitating rapid and robust amplification of short CyHV-3-TK amplicons. This strategic design not only enhances assay sensitivity but also accelerates reaction kinetics, as shorter amplicons enable faster and more efficient amplification (Daher et al. 2016). Moreover, the specificity of the crRNA is a critical determinant of target recognition kinetics. In this study, the CyHV-3-specific crRNA sequence was designed based on a conserved region of the CyHV-3-TK gene within the RAA-amplified product. Sequence overlap with RAA primers was avoided to minimise the risk of false-positive signals (Kellner et al. 2019). The SHERLOCK assay leverages LwaCas13a's RNA-guided RNA-targeting capability, wherein the enzyme, upon activation by target recognition, indiscriminately cleaves nearby RNA molecules, enabling signal amplification through reporter degradation (East-Seletsky et al. 2016). Particularly, LwaCas13a does not exhibit stringent sequence constraints at the target site, allowing for greater flexibility in target selection and expanding the range of detectable sequences (Gootenberg et al. 2018).

The results of this study demonstrated that the KHV-SHERLOCK assay exhibits significantly higher sensitivity than conventional PCR targeting the same CyHV-3-TK gene. Specifically, the KHV-SHERLOCK assay detected as little as 21 CyHV-3-TK plasmid copies, whereas conventional PCR required 210 copies. This enhanced sensitivity can be attributed to the integration of RAA, which exponentially amplifies trace amounts of target DNA, thereby increasing the substrate availability for LwaCas13a-mediated detection. The findings are consistent with previous studies reporting the ultra-sensitivity (2 aM to 8 zM) of SHERLOCK-based assays (Gootenberg et al. 2018; Kellner et al. 2019; Zhao et al. 2022; Li, Zhang, et al. 2024; Li, Li, et al. 2024; Wang et al. 2025). The combination of isothermal RAA preamplification and CRISPR-Cas13a collateral cleavage activity ensures robust detection of low-abundance targets,

establishing SHERLOCK as a powerful diagnostic tool for CyHV-3 surveillance.

The KHV-SHERLOCK assay exhibited complete concordance with conventional PCR-positive samples, confirming its high diagnostic sensitivity and supporting its utility in clinical settings for reliable detection of CyHV-3.

In conclusion, the integrated assay combining RAA with a CRISPR-Cas13a system and LFD detection enables rapid, visual identification of CyHV-3, offering significant advantages over conventional molecular diagnostic methods for POC applications. This protocol demonstrates operational simplicity, clear visual interpretation of results, and high specificity and sensitivity, making it a promising tool for on-site diagnosis of CyHV-3 infections in clinical settings.

Author Contributions

H.S. and M.S. conceived and designed the study. H.S. performed the experiments. N.A. helped with the experiments. H.S. analyzed the data and wrote the manuscript. M.S. and N.A. revised the manuscript. All authors read and approved the final manuscript.

Acknowledgements

The Hatem Soliman gratefully acknowledges the financial support provided by the Academy of Scientific Research and Technology (ASRT), Egypt, through the ASRT-STAR fellowship program, Cycle 1. This work was supported in part by the Veterinary Medicine University, Vienna. Open Access funding provided by Veterinarmedizinische Universität Wien/KEMÖ.

Conflicts of Interest

The authors declare no conflicts of interest.

Data Availability Statement

Data sharing not applicable to this article as no datasets were generated or analysed during the current study.

References

Abudayyeh, O. O., J. S. Gootenberg, S. Konermann, et al. 2016. "C2c2 Is a Single-Component Programmable RNA-Guided RNA-Targeting CRISPR Effector." *Science* 353, no. 6299: aaf5573.

Allan-Blitz, L.-T., P. Shah, G. Adams, et al. 2023. "Development of Cas13a-Based Assays for *Neisseria gonorrhoeae* Detection and Gyrase A Determination." *mSphere* 8, no. 5: e0041623.

Aoki, T., I. Hirono, K. Kurokawa, et al. 2007. "Genome Sequences of Three Koi Herpesvirus Isolates Representing the Expanding Distribution of an Emerging Disease Threatening Koi and Common Carp Worldwide." *Journal of Virology* 81: 5058–5065.

Bavarsad, M., A. Abed-Elmdoost, M. R. Tabandeh, M. Alishahi, A. Mirvaghefi, and H. Hamid Farahmand. 2025. "Rapid Cyprinid Herpesvirus 3 Detection Using Loop-Mediated Isothermal Amplification (LAMP) Combined With Gold Nanoprobes and SYBR Safe, Aquaculture, 594, 741374." *Journal of Aquaculture* 2024: 741374.

Bei, L., C. Hairong, Y. Qingfeng, et al. 2010. "Rapid Amplification of Nucleic Acids Using Recombinase-Mediated Amplification Technology." *Scientia Sinica (Vita)* 40, no. 10: 983–988.

Bercovier, H., Y. Fishman, R. Nahary, et al. 2005. "Cloning of the Koi Herpesvirus (KHV) Gene Encoding Thymidine Kinase and Its Use for a Highly Sensitive PCR Based Diagnosis." *BMC Microbiology* 5: 1–9.

Bergmann, S. M., M. Riechardt, D. Fichtner, P. Lee, and J. Kempter. 2010. "Investigation on the Diagnostic Sensitivity of Molecular Tools Used for Detection of Koi Herpesvirus." *Journal of Virological Methods* 163: 229–233.

Boivin, G., Z. Coulombe, and S. Rivest. 2002. "Intranasal Herpes Simplex Virus Type 2 Inoculation Causes a Profound Thymidine Kinase Dependent Cerebral Inflammatory Response in the Mouse Hindbrain." *European Journal of Neuroscience* 16: 29–43.

Calderón, I. L., M. J. Barros, N. Fernández-Navarro, and L. G. Acuña. 2024. "Detection of Nucleic Acids of the Fish Pathogen *Yersinia ruckeri* From Planktonic and Biofilm Samples With a CRISPR/Cas13a-Based Assay." *Microorganisms* 12, no. 2: 283.

Chen, J. S., E. Ma, L. B. Harrington, et al. 2018. "CRISPR-Cas12a Target Binding Unleashes Indiscriminate Single-Stranded DNase Activity." *Science* 360: 436–439.

Daher, R. K., G. Stewart, M. Boissinot, and M. G. Bergeron. 2016. "Recombinase Polymerase Amplification for Diagnostic Applications." *Clinical Chemistry* 62: 947–958.

East-Seletsky, A., M. R. O'Connell, S. C. Knight, et al. 2016. "Two Distinct RNase Activities of CRISPR-C2c2 Enable Guide-RNA Processing and RNA Detection." *Nature* 538: 270–273.

El-Matbouli, M., U. Rucker, and H. Soliman. 2007. "Detection of Cyprinid Herpesvirus-3 (CyHV-3) DNA in Infected Fish Tissues by Nested Polymerase Chain Reaction." *Diseases of Aquatic Organisms* 78: 23–28.

Engelsma, M. Y., K. Way, M. J. Dodge, et al. 2013. "Detection of Novel Strains of Cyprinid Herpesvirus Closely Related to Koi Herpesvirus." *Diseases of Aquatic Organisms* 107: 113–120.

Fan, G. H., X. X. Shen, F. Li, et al. 2019. "Development of an Internally Controlled Reverse Transcription Recombinase-Aided Amplification Assay for the Rapid and Visual Detection of West Nile Virus." *Biomedical and Environmental Sciences* 32: 926–929.

Gilad, O., S. Yun, M. A. Adkison, et al. 2003. "Molecular Comparison of Isolates of an Emerging Fish Pathogen, Koi Herpesvirus, and the Effect of Water Temperature on Mortality of Experimentally Infected Koi." *Journal of General Virology* 84: 2661–2668.

Gilad, O., S. Yun, F. J. Zagmutt-Vergara, C. M. Leutenegger, H. Bercovier, and R. P. Hedrick. 2004. "Concentrations of a Koi Herpesvirus (KHV) in Tissues of Experimentally Infected *Cyprinus carpio koi* as Assessed by Real-Time TaqMan PCR." *Diseases of Aquatic Organisms* 60: 179–187.

Gootenberg, J. S., O. O. Abudayyeh, M. J. Kellner, J. Joung, J. J. Collins, and F. Zhang. 2018. "Multiplexed and Portable Nucleic Acid Detection Platform With Cas13, Cas12a, and Csm6." *Science* 360: 439–444.

Gootenberg, J. S., O. O. Abudayyeh, J. W. Lee, et al. 2017. "Nucleic Acid Detection With CRISPR-Cas13a/C2c2." *Science* 356: 438–442.

Gotesman, M., J. Kattlun, S. M. Bergmann, and M. El-Matbouli. 2013. "CyHV-3: The Third Cyprinid Herpesvirus." *Diseases of Aquatic Organisms* 105, no. 2: 163–174. <https://doi.org/10.3354/dao02614>.

Gui, L., Y. Zhao, D. Xu, et al. 2022. "Quick Detection of *Carassius auratus* Herpesvirus (CaHV) by Recombinase-Aid Amplification Lateral Flow Dipstick (RAA-LFD) Method." *Frontiers in Cellular and Infection Microbiology* 12: 981911.

Gunimaladevi, I., T. Kono, M. Venugopal, and M. Sakai. 2004. "Detection of Koi Herpesvirus in Common Carp, *Cyprinus carpio* L., by Loop-Mediated Isothermal Amplification." *Journal of Fish Diseases* 27: 583–589.

- Haenen, O., K. Way, S. Bergmann, and E. Ariel. 2004. "The Emergence of Koi Herpesvirus and Its Significance to European Aquaculture." *Bulletin of the European Association of Fish Pathology* 24: 293–307.
- Hatoum-Aslan, A. 2018. "CRISPR Methods for Nucleic Acid Detection Herald the Future of Molecular Diagnostics." *Clinical Chemistry* 64: 1681–1683.
- Hedrick, R. P., O. Gilad, S. Yun, et al. 2000. "A Herpesvirus Associated With Mass Mortality of Juvenile and Adult Koi, a Strain of Common Carp." *Journal of Aquatic Animal Health* 12: 44–57.
- ICTV. 2022. "International Committee on Taxonomy of Viruses—Virus Taxonomy Release." <https://ictv.global/taxonomy>.
- Kellner, M. J., J. G. Koob, J. S. Gootenberg, O. O. Abudayyeh, and F. Zhang. 2019. "SHERLOCK: Nucleic Acid Detection With CRISPR Nucleases." *Nature Protocols* 14, no. 10: 2986–3012.
- Lee, R. A., H. Puig, P. Q. Nguyen, et al. 2020. "Ultrasensitive CRISPR-Based Diagnostic for Field-Applicable Detection of *Plasmodium* Species in Symptomatic and Asymptomatic Malaria." *Proceedings of the National Academy of Sciences United States of America* 117, no. 41: 25722–25731.
- Li, J., Z. Zhang, Z. Zhang, et al. 2024. "Rapid Detection of Avian Leukemia Virus Using CRISPR/Cas13a Based Lateral Flow Dipstick." *Frontiers in Veterinary Science* 11: 1424238.
- Li, L., S. Li, N. Wu, et al. 2019. "HOLMESv2: A CRISPR- Cas12b-Assisted Platform for Nucleic Acid Detection and DNA Methylation Quantitation." *ACS Synthetic Biology Journal* 8: 2228–2237.
- Li, S. Y., Q. X. Cheng, J. K. Liu, X. Q. Nie, G. P. Zhao, and J. Wang. 2018. "CRISPR-Cas12a Has Both Cis- and Trans-Cleavage Activities on Single-Stranded DNA." *Cell Research* 28: 491–493.
- Li, X., S. Zhu, X. Zhang, et al. 2023. "Advances in the Application of Recombinase-Aided Amplification Combined With CRISPR-Cas Technology in Quick Detection of Pathogenic Microbes." *Frontiers in Bioengineering and Biotechnology* 11: 1215466.
- Li, Y., R. Li, X. Mo, et al. 2024. "Development of Real-Time Recombinase Polymerase Amplification (RPA) and RPA Combined With Lateral Flow Dipstick (LFD) Assays for the Rapid and Sensitive Detection of Cyprinid Herpesvirus 3." *Journal of Fish Diseases* 47: e13960.
- Myhrvold, C., C. A. Freije, J. S. Gootenberg, et al. 2018. "Field-Deployable Viral Diagnostics Using CRISPR-Cas13." *Science* 360: 444–448.
- Perelberg, A., M. Smirnov, M. Hutoran, A. Diamant, Y. Bejerano, and M. Kotler. 2003. "Epidemiological Description of a New Viral Disease Afflicting Cultured *Cyprinus carpio* in Israel." *Israeli Journal of Aquaculture* 55: 5–12.
- Prescott, M. A., A. N. Reed, L. Jin, and M. K. Pastey. 2016. "Rapid Detection of Cyprinid Herpesvirus 3 in Latently Infected Koi by Recombinase Polymerase Amplification." *Journal of Aquatic Animal Health* 28: 173–180.
- Ramadan, N. K., N. Gaber, N. M. Ali, O. S. O. Amer, and H. Soliman. 2025. "SHERLOCK, a Novel CRISPR-Cas13a-Based Assay for Detection of Infectious Bursal Disease Virus." *Journal of Virological Methods* 337: 115185.
- Ronen, A., A. Perelberg, J. Abramowitz, et al. 2003. "Efficient Vaccine Against the Virus Causing a Lethal Disease in Cultured *Cyprinus carpio*." *Vaccine* 21: 4677–4684.
- Saleh, M., and M. El-Matbouli. 2015. "Rapid Detection of Cyprinid Herpesvirus-3 (CyHV-3) Using a Gold Nanoparticle-Based Hybridization Assay." *Journal of Virological Methods* 217: 50–54. <https://doi.org/10.1016/j.jviromet.2015.02.021>.
- Soliman, H., and M. El-Matbouli. 2005. "An Inexpensive and Rapid Diagnostic Method of Koi Herpesvirus (KHV) Infection by Loop-Mediated Isothermal Amplification." *Virology Journal* 2: 1–8.
- Soliman, H., and M. El-Matbouli. 2010. "Loop Mediated Isothermal Amplification Combined With Nucleic Acid Lateral Flow Strip for Diagnosis of Cyprinid Herpesvirus-3." *Molecular and Cellular Probes* 24: 38–43.
- Soliman, H., and M. El-Matbouli. 2018. "Rapid Detection and Differentiation of Carp Oedema Virus and Cyprinid Herpes Virus-3 in Koi and Common Carp." *Journal of Fish Diseases* 41: 761–772.
- Sullivan, T. J., A. K. Dhar, R. Cruz-Flores, and A. G. Bodnar. 2019. "Rapid, CRISPR-Based, Field-Deployable Detection of White Spot Syndrome Virus in Shrimp." *Scientific Reports* 9: 19702.
- Thompson, J. D., D. G. Higgins, and T. J. Gibson. 1994. "CLUSTAL W: Improving the Sensitivity of Progressive Multiple Sequence Alignment Through Sequence Weighting, Position-Specific Gap Penalties and Weight Matrix Choice." *Nucleic Acids Research* 22: 4673–4680.
- Waltzek, T. B., G. O. Kelley, D. M. Stone, et al. 2005. "Koi Herpesvirus Represents a Third Cyprinid Herpesvirus (CyHV-3) in the Family Herpesviridae." *Journal of General Virology* 86, no. Pt 6: 1659–1667. <https://doi.org/10.1099/vir.0.80982-0>.
- Wang, H., Y. Yan, C. R. Codjia, et al. 2025. "CRISPR-Cas13a Coupled With Recombinase Aided Amplification: Development of an Ultrasensitive Point-of-Care Detection Method for Visual Diagnosis of Egg Drop Syndrome Virus." *Talanta* 296: 128470.
- Wang, K., R. Huang, L. Zhang, D. Liu, and Y. Diao. 2022. "Recombinase-Aided Amplification Combined With Lateral Flow (LF-RAA) Assay for Rapid AAV Genome Detection." *ACS Omega Journal* 7, no. 51: 47832–47839.
- Wang, R. H., H. Zhang, Y. Zhang, et al. 2019. "Development and Evaluation of Recombinase-Aided Amplification Assays Incorporating Competitive Internal Controls for Detection of Human Adenovirus Serotypes 3 and 7." *Virology Journal* 16: 86.
- World Organization for Animal Health. 2022. *Manual of Diagnostic Tests for Aquatic Animals. Chapter 2.3.6. Infection with Koi Herpesvirus*. WOAH. <https://www.woah.org/en/disease/koi-herpesvirus-disease/>. (Accessed 15.06.2025).
- Yan, C., Y. Zhou, S. Du, et al. 2023. "Recombinase-Aided Amplification Assay for Rapid Detection of Hypervirulent *Klebsiella pneumoniae* (hvKp) and Characterization of the hvKp Pathotype." *Microbiology Spectrum Journal* 11, no. 2: e0398422.
- Yin, L. J., S. L. Man, S. Y. Ye, G. Z. Liu, and L. Ma. 2021. "CRISPR-Cas Based Virus Detection: Recent Advances and Perspectives." *Biosensors and Bioelectronics* 193: 113541.
- Yoshino, M., H. Watari, T. Kojima, M. Ikeda, and J. Kurita. 2009. "Rapid, Sensitive and Simple Detection Method for Koi Herpesvirus Using Loop-Mediated Isothermal Amplification." *Microbiology and Immunology* 53: 375–383.
- Yuasa, K., M. Sano, J. Kurita, T. Ito, and T. Iida. 2005. "Improvement of a PCR Method With the Sph1-5 Primer Set for the Detection of Koi Herpesvirus (KHV)." *Fish Pathology* 40: 37–39.
- Zhao, J., Y. Li, Q. Xue, Z. Zhu, M. Zou, and F. Fang. 2022. "A Novel Rapid Visual Detection Assay for *Toxoplasma Gondii* Combining Recombinase-Aided Amplification and Lateral Flow Dipstick Coupled With CRISPR-Cas13a Fluorescence (RAA-Cas13a-LFD)." *Parasite* 29: 21.
- Zhao, S., Q. Zhang, X. Wang, et al. 2023. "Development and Performance of Recombinase-Aided Amplification (RAA) Assay for Detecting *Schistosoma haematobium* DNA in Urine Samples." *Heliyon* 9, no. 12: e23031.